(en)

Composition
easy-graft CLASSIC/CLASSIC* contains:
 • Pure beta-tricalcium phosphate (β-TCP) granules coated with poly(lactide-co-glycolide) (PLGA)
 • BioLinker (N-Methyl-2-pyrrolidone solution)
Description
easy-graft CLASSIC/CLASSIC* is a sterile, synthetic bone defect filler. It consists of two components: granules (syringe) and BioLinker (ampulla). After mixing, *easy-graft* CLASSIC/CLASSIC* is putty-like and can be applied directly from the syringe into the bone defect. In contact with blood or other body fluids, *easy-graft* CLASSIC/CLASSIC* hardens within minutes to form a solid but porous implant.
easy-graft CLASSIC/CLASSIC* is biocompatible and osteo-regenerative. For β-TCP and PLGA in adequate form, no inflammatory or immunoreactions are known.
easy-graft CLASSIC/CLASSIC* forms a stable and porous body within the defect. The granules already feature a high porosity. After placing the material in the defect, the pore volume allows the uptake of blood. With the integration of blood, important factors for tissue regeneration are incorporated. In a second step, autogenous tissue can grow into the pore structure.
 During the initial phase of degradation, *easy-graft* CLASSIC/CLASSIC* may swell by taking up body fluids, thus supporting a close contact to the surrounding bone tissue. When applied in larger defects, this might result in slight sensation of pressure by the patient.
easy-graft CLASSIC/CLASSIC* contains no animal or human derived substances. β-TCP and PLGA are derived from fully synthetic raw materials.

Indications
easy-graft CLASSIC/CLASSIC* should be used for the filling of mostly unloaded and preferably multi-walled dental or maxillo-facial bone defects. The implantation site should be free of infections and soft or granulation tissue. Indications may include:
 • Defects after removal of bone cysts
 • Periodontal defects
 • Peri-implant defects
 • Augmentation of alveolar crest (e.g. Guided Bone Regeneration, GBR)
 • Defects after root end resection
 • Extraction defects
 • Defects after surgical removal of retained teeth
 • Sinus floor elevations
 • Defects after removal of autologous bone

Application
 Aseptic working methods are essential during application.
easy-graft CLASSIC are spherical shaped granules provided in particle sizes of 500-630 µm for smaller defects and 500-1000 µm for larger defects.
easy-graft CLASSIC* are irregularly shaped granules provided in the particle sizes of 450-630 µm for smaller defects and 450-1000 µm for larger defects.
 The amount of *easy-graft* CLASSIC/CLASSIC* used per patient during one appointment or surgical application should be restricted to a total of 1.8 ml.
 For the application of *easy-graft* CLASSIC/CLASSIC* refer to the instructions for use table (steps I. to VI.).
 Dental membranes may be used in combination with *easy-graft* CLASSIC/CLASSIC* at the discretion of the practitioner.

Contraindications
easy-graft CLASSIC/CLASSIC* must not be used during pregnancy.
easy-graft CLASSIC/CLASSIC* should not be applied in case of the presence of any diseases or therapies which are detrimental to the healing of bone.
 No allergic reactions against *easy-graft* CLASSIC/CLASSIC* or its components are known. However, if an allergic reaction is suspected, the product should not be applied in the affected patient.

Precautions
 In order to facilitate the formation of new bone to a maximum extent, *easy-graft* CLASSIC/CLASSIC* should only be implanted in close contact to well vascularized surrounding bony tissue (refreshment of bleeding might be necessary).
 Before filling periodontal defects, a sufficient pre-treatment is necessary. This should include sufficient plaque control and a successful local treatment of the periodontal lesion (e.g. root planning, debridement).

Possible adverse effects
 In general, possible adverse reactions are comparable to those encountered in autogenous bone or bone graft substitute grafting procedures. These include, but are not limited to: edema, bleeding after surgery, numbness by injured nerve, infection of the soft tissue and/or bone, soft tissue ulceration, incomplete bone formation, lack of bone formation, fracture or inflammation of the newly formed bone.
 Allergic or inflammatory reactions to the implant material or its decomposition products have not been recorded.

Interaction with other substances
 Not known.


Warnings
easy-graft CLASSIC/CLASSIC* is restricted to dental and maxillo-facial applications and to be sold to, or on the order of a licensed physician or dentist.
easy-graft CLASSIC/CLASSIC* must not be applied to acute infected defects. The implantation site should be free of infections before applying the *easy-graft* CLASSIC/CLASSIC*.
 Do not place *easy-graft* CLASSIC/CLASSIC* in direct contact with nerve tissue. The product has not been clinically tested in children.
 Lactation: It is not known whether metabolites of *easy-graft* CLASSIC/CLASSIC* pass into breast milk.

High N-methyl-2-pyrrolidone (contained in the BioLinker activator) exposure may cause malaise, headache, stomach pain and nausea. Therefore, the amount of *easy-graft* CLASSIC/CLASSIC* used per patient during one appointment or surgical application should be restricted to a total of 1.8 ml.
 Do not use if package is opened or damaged, or if the expiration date has been exceeded.

easy-graft CLASSIC/CLASSIC* cannot be reused. Discard any unused material. Sterility and function are no longer guaranteed in the event of re-use.
Properties
easy-graft CLASSIC/CLASSIC* consists of β-TCP, PLGA and BioLinker. Resorption proceeds via three steps: 1) The BioLinker is metabolized and excreted within days mainly through the urine. 2) The absorbable polymer (PLGA) is cleaved by hydrolysis. The degradation products are excreted. 3) The resorption of β-TCP takes place mostly parallel to bone regeneration. Depending on the regeneration potential of the tissue, *easy-graft* CLASSIC/CLASSIC* will be completely resorbed within 5 to 15 months.
easy-graft CLASSIC/CLASSIC* supports bone regeneration and bone formation. The application of the material does not guarantee a treatment success. The success of the therapy depends on many factors such as surgical technique and habits, age and the regeneration potential of the patient.

Miscellaneous
easy-graft CLASSIC/CLASSIC* is gamma sterilized and intended for single application. Unused material can be disposed in the domestic waste.
Storage instructions
easy-graft CLASSIC/CLASSIC* must be stored in its original packaging. Store between 5°C/41°F to 25°C/77°F. Keep away from sunlight.
Manufacturer
 Degradable Solutions AG, Wagistrasse 23, 8952 Schlieren / Zurich, Switzerland

I.


Open ampulla with BioLinker by twisting off the lid, subsequently open the syringe with granules by pulling out the plunger at the rear of the syringe.
Caution: Plunger will be needed again, omit contact with contaminated surfaces.
 By holding syringe with the opened end upwards, the loss of granules can be avoided.

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II.


Add BioLinker through the funnel-shaped opening into the syringe to the granules. Pulling back the plug of the syringe slightly facilitates the wetting of the granules.

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III.

Place the plunger back into the syringe at the rear end (a). The granules must be completely wetted with BioLinker. In case of incomplete wetting of the granules, a complete wetting can be achieved by moving back and forth the plunger and the plug 1-3 times.
 Then the plug is removed (b).

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IV.


Before application, remove spare BioLinker by pushing *easy-graft* CLASSIC/CLASSIC* to the very front of syringe. The mass should not be compressed during this step, in order not to compromise smooth ejection of the material from the syringe.

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
VI.

The material may be shaped and compressed in the defect. In contact with blood or other body fluids, *easy-graft* CLASSIC/CLASSIC* hardens within minutes to form a solid but porous implant.

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Degradable Solutions AG, Wagistrasse 23, 8952 Schlieren / Zurich, Suisse