GUIDOR® *easy-graft*® CRYSTAL⁺ en

in situ hardening synthetic bone graft system

Composition
easy-graft CRYSTAL⁺ contains:

- Biphasic calcium phosphate (60 % hydroxyapatite, 40 % beta-tricalcium phosphate) granules coated with poly(lactide-co-glycolide) (PLGA)
- BioLinker (N-Methyl-2-pyrrolidone solution)

Description
easy-graft CRYSTAL⁺ is a sterile, synthetic bone defect filler. It consists of two components: granules (syringe) and BioLinker (ampulla). After mixing, *easy-graft* CRYSTAL⁺ is putty-like and can be applied directly from the syringe into the bone defect. In contact with blood or other body fluids, *easy-graft* CRYSTAL⁺ hardens within minutes to form a solid but porous implant.

easy-graft CRYSTAL⁺ is biocompatible and osteo-regenerative. For biphasic calcium phosphate and PLGA in adequate form, no inflammatory or immunoreactions are known. *easy-graft* CRYSTAL⁺ forms a stable and porous body within the defect. The granules feature a high porosity. After placing the material in the defect, the pore volume allows the uptake of blood. With the integration of blood, important factors for tissue regeneration are incorporated. In a second step, autogenous tissue can grow into the pore structure. During the initial phase of degradation, *easy-graft* CRYSTAL⁺ may swell by taking up body fluids, thus supporting a close contact to the surrounding bone tissue. When applied in larger defects, this might result in slight sensation of pressure by the patient. *easy-graft* CRYSTAL⁺ contains no animal- or human-derived substances. Biphasic calcium phosphate and PLGA are derived from fully synthetic raw materials.

Indications
easy-graft CRYSTAL⁺ should be used for the filling of mostly unloaded and preferably multi-walled dental or maxillo-facial bone defects. The implantation site should be free of infections and soft or granulation tissue. Indications may include:

- Defects after removal of bone cysts
- Periodontal defects
- Peri-implant defects
- Augmentation of alveolar crest (e.g. Guided Bone Regeneration, GBR)
- Defects after root end resection
- Extraction defects
- Defects after surgical removal of retained teeth
- Sinus floor elevations
- Defects after removal of autologous bone

Application
Aseptic working methods are essential during application. *easy-graft* CRYSTAL⁺ are irregularly shaped granules provided in the particle sizes of 450 – 630 µm for smaller defects and 450 – 1000 µm for larger defects. The amount of *easy-graft* CRYSTAL⁺ used per patient during one appointment or surgical application should be restricted to a total of 1.8 ml. For the application of *easy-graft* CRYSTAL⁺ refer to the instructions for use table (steps I. to VI.). Dental membranes may be used in combination with *easy-graft* CRYSTAL⁺ at the discretion of the practitioner.

Contraindications
easy-graft CRYSTAL⁺ must not be used during pregnancy. *easy-graft* CRYSTAL⁺ should not be applied in case of the presence of any diseases or therapies which are detrimental to the healing of bone. No allergic reactions against *easy-graft* CRYSTAL⁺ or its components are known. However, if an allergic reaction is suspected, the product should not be applied in the affected patient.

Precautions
In order to facilitate the formation of new bone to a maximum extent, *easy-graft* CRYSTAL⁺ should only be implanted in close contact to well-vascularized surrounding bony tissue (re-freshment of bleeding might be necessary). Care should be taken that *easy-graft* CRYSTAL⁺ granules are not introduced into surrounding soft tissue, because the granules will only partially resorb. Before filling periodontal defects, a sufficient pre-treatment is necessary. This should include sufficient plaque control and a successful local treatment of the periodontal lesion (e.g. root planing, debridement).

Possible adverse effects
In general, possible adverse reactions are comparable to those encountered in autogenous bone or bone graft substitute grafting procedures. These include, but are not limited to: edema, bleeding after surgery, numbness by injured nerve, infection of the soft tissue and/or bone, soft tissue ulceration, incomplete bone formation, lack of bone formation, fracture or inflammation of the newly formed bone. Allergic or inflammatory reactions to the implant material or its decomposition products have not been recorded.

Interaction with other substances
Not known.

Warnings
easy-graft CRYSTAL⁺ is restricted to dental and maxillo-facial applications and to be sold to, or on the order of a licensed physician or dentist. *easy-graft* CRYSTAL⁺ must not be applied to acute infected defects. The implantation site should be free of infections before applying the *easy-graft* CRYSTAL⁺. Do not place *easy-graft* CRYSTAL⁺ in direct contact with nerve tissue. The product has not been clinically tested in children. Lactation: It is not known whether metabolites of *easy-graft* CRYSTAL⁺ pass into breast milk.

BioLinker contains N-Methyl-2-pyrrolidone. N-Methyl-2-pyrrolidone is irritant to eyes, respiratory system and skin. Wear gloves and protection goggles when working with BioLinker. Change gloves after contact with BioLinker. After eye contact: In case of contact, immediately flush eyes with copious amounts of water for at least 10 minutes and consult a doctor.

High N-methyl-2-pyrrolidone (contained in the BioLinker activator) exposure may cause malaise, headache, stomach pain and nausea. Therefore, the amount of *easy-graft* CRYSTAL⁺ used per patient during one appointment or surgical application should be restricted to a total of 1.8 ml.

Do not use if package is opened or damaged, or if the expiration date has been exceeded. *easy-graft* CRYSTAL⁺ cannot be reused. Discard any unused material. Sterility and function are no longer guaranteed in the event of re-use.

Properties
easy-graft CRYSTAL⁺ consists of biphasic calcium phosphate, PLGA and BioLinker. Resorption proceeds via three steps: 1) The BioLinker is metabolized and excreted within days mainly through the urine. 2) The absorbable polymer (PLGA) is cleaved by hydrolysis. The degradation products are excreted. 3) The biphasic calcium phosphate is partially resorbable by cellular mechanisms and direct dissolution. A part of the biphasic calcium phosphate remains embedded in autogenous tissue for prolonged periods of time in order to preserve the volume. *easy-graft* CRYSTAL⁺ supports bone regeneration and bone formation. The application of the material does not guarantee a treatment success. The success of the therapy depends on many factors such as surgical technique and habits, age and the regeneration potential of the patient.

Miscellaneous
easy-graft CRYSTAL⁺ is gamma sterilized and intended for single application. Unused material can be disposed in the domestic waste.

Storage instructions
easy-graft CRYSTAL⁺ must be stored in its original packaging. Store between 5°C/41°F to 25°C/77°F. Keep away from sunlight.


Manufacturer
Degradable Solutions AG, Wagistrasse 23, 8952 Schlieren / Zurich, Switzerland

CE 0297	CE Marking - Notified body: DQS Medizinprodukte GmbH
REF	Catalogue number
ⓧ	Do not re-use
LOT	Batch code
ⓘ	Use-by date
ⓘ	Consult instructions for use
STERILIS R	Sterilized using irradiation
15°C - 25°C	Temperature limit
☀	Keep away from sunlight
☹	Do not use if package is damaged
🏭	Manufacturer
MD	Medical Device
EC/REP	European Authorized Representative

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I.


Open ampulla with BioLinker by twisting off the lid, subsequently open the syringe with granules by pulling out the plunger at the rear of the syringe.
Caution: Plunger will be needed again, omit contact with contaminated surfaces.
 By holding syringe with the opened end upwards, the loss of granules can be avoided.

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II.


Add BioLinker through the funnel-shaped opening into the syringe to the granules. Pulling back the plug of the syringe slightly facilitates the wetting of the granules. *easy-graft* CRYSTAL⁺ should be applied within one minute after the addition of BioLinker.

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III.


Place the plunger back into the syringe at the rear end (a). The granules must be completely wetted with BioLinker. In case of incomplete wetting of the granules, a complete wetting can be achieved by moving back and forth the plunger and the plug 1-3 times. Then the plug is removed (b).

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IV.


Before application, remove spare BioLinker by pushing *easy-graft* CRYSTAL⁺ to the very front of syringe. The mass should not be compressed during this step in order not to compromise moist ejection of the material from the syringe.

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V.

Apply *easy-graft* CRYSTAL⁺ directly from the syringe into the prepared defect. **Caution:** The amount of *easy-graft* CRYSTAL⁺ used per patient during one appointment or surgical application should be restricted to a total of 1.8 ml.

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VI.

The material may be shaped and compressed in the defect. In contact with blood or other body fluids, *easy-graft* CRYSTAL⁺ hardens within minutes to form a solid but porous implant.

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